

To: Zimmer Dental Customers


Date: June 23, 2006

From: Zimmer Dental Marketing Department


Subject: Modification to the Healing Collars for *Tapered Screw-Vent*® Implants

Zimmer Dental has made a minor modification to the Healing Collars for *Tapered Screw-Vent* implants.

The wide section of the shaft above the screw threads has been narrowed on 3.5mm sizes and shortened on 4.5mm and 5.7mm sizes. The illustration below shows a before (rev T) and after (rev V) image of a 3.5mm Healing Collar. This part of the healing collar does not come in contact with the implant and this change has no affect on the function of the part.



REV 'V'



REV 'T'

If you have any questions, please contact our Technical Services Department at 1+760-929-4300.

Thank you for choosing Zimmer Dental.