




Product Information

New 4.1mmD Tapered Screw-Vent® Implant Dimensions




4.1mmD Tapered Screw-Vent Drilling Sequence


4.1mmD Tapered Screw-Vent Implant (3.5mmD Platform)




1
SV2.3D
2.3mmD
Drill




2
SV2.8D
2.8mmD
Drill



3
FOR SOFT BONE
SV3.4D
3.4mmD
Drill



3
FOR DENSE BONE
TSV3.8D
3.8mm/3.4mmD
Drill



4
OPTIONAL FOR DENSE BONE
TT4.1
4.1mmD
Cortical Bone Tap